

Evaluation of 3-D grain motion inside a granular assembly during shear deformation

T. Matsushima (9483)¹⁾, A. Tsuchiyama (4116)²⁾, K. Uesugi (1544)³⁾, T. Nakano (4129)⁴⁾, K. Toda (9495)⁵⁾, Y. Kawamura (9487)⁵⁾, K. Sato (14013)⁶⁾

(1) Institute of Engineering Mechanics and Systems, University of Tsukuba, (2) Department of Earth and Space Science, Osaka University, (3) Spring-8/JASRI, (4) Geological Survey of Japan/AIST and (5) Graduate School of Engineering, University of Tsukuba. (6) College of Engineering Systems, University of Tsukuba

Mechanics of granular materials is used in various engineering fields such as geotechnical engineering and powder engineering. Within a quasi-static regime, a granular assembly stands against an external force by constructing a contact-force network. Therefore it is essential to study its packing structure (number of contact points, their orientation etc.) and its evolution during shear deformation. This research attempted such particle-level observation inside a triaxial specimen using micro X-ray CT system (SP- μ CT) at BL20B2.

We developed a micro triaxial apparatus suitable for SP- μ CT, which enables us to keep a confining pressure constant and to apply a precise axial strain to a small specimen with 3.4 mm in diameter and 10 mm in height. The material used is Toyoura sand, a standard sand commonly used for geotechnical engineering research, whose grain size ranges from 0.1 to 0.2 mm. Two types of specimen, a loose one and a dense one were tested and the stress - strain curves were obtained (Fig.1). The

resultant peak axial stresses for both specimens are rather higher than those with a conventional tri-axial test whose specimen size is over 50 mm in diameter and 100mm in height. It is mainly due to the effect of relatively-thick (0.15mm) membrane wrapping the specimens.

Loading was stopped several times during the test, as it is observed in Fig.1 by the vertical stress drop, to take a series of X-ray images. Fig.2 shows the reconstructed CT images of a central vertical cross section of the dense specimen before and after loading. It can be recognized that the shear deformation is concentrated on the upper part of the specimen, accompanying a considerable dilation (volume increase). The images are clear enough to quantify each grain motion including its rotation. Such information as 3-D microscopic motion of grains in a shear zone, which has never been observed before, is quite useful for the verification of microscopic constitutive model as well as particle-level simulation.




Fig.1 Stress-strain curve




Fig.2 Reconstructed CT images.
 $\varepsilon_a = 0$ (left) and $\varepsilon_a = 15\%$ (right)